GE Healthcare

Revolution™ EVO

User Manual

GE Healthcare Japan Corporation does business as GE Healthcare

This manual supports the following configurations:

- Revolution™ EVO EL
- Revolution™ EVO EX
- Revolution™ EVO ES
- Revolution™ EVO
- Revolution™ EVO 48kW

Not all configurations are available in all regions.

This product is certified as a Revolution™ EVO CT Scanner.

Revolution™ EVO
Operator Manual, English
5805440-1EN
Revision: 1
© 2019 General Electric Company
All rights reserved.

SECRET DE AFACERI

CONFIDENTIAL

231

Gantry

Gantry display

The gantry display provides gantry and table status information.


Figure 7-9: Gantry display with patient information

Table 7-7: Gantry display components

No.	Component	Description
1.	Collision indicator 	The Collision light illuminates and a picture of CT system with high-lighted collision sensors is posted on the gantry display (Figure 7-10) when there is a possibility that the table, cradle and gantry comes in contact with each other. There are tilt/elevation collision sensors on the front and rear covers at the top of the gantry opening. The collision sensors are active during table elevation and during tilting gantry. The light also illuminates if there is 4.5 kg (10 lb) of resistance during cradle motion in/out or when reaching the travel limits of table, cradle and gantry motion. In the case of interference, you may need to raise or lower the table, check the gantry tilt, or determine if resistance is being caused by any patient restraints, accessories, or by the patient. If you need to tilt the gantry after a collision has been detected, the gantry only tilts in 0.5 degree increments.  The picture of CT system with high-lighted collision sensors on the gantry display can be hidden by touching the display.
2.	Cradle unlocked indicator 	The Cradle Unlocked Indicator is illuminated, indicating the cradle is unlocked. When the cradle is unlocked, it is "free floating".
3.	Laser light indicator	The Laser Light Indicator is illuminated when the laser alignment

CONFIDENTIAL

No.	Component	Description
		light is on. Normally, this light is only on when positioning a patient.  CAUTION: For patient safety, it is important to always have patients close their eyes anytime the laser alignment light is on.
4.	Lung indicator 	The lung indicator is not lit on this system.
5.	Heart indicator 	The heart is lit whenever there is an ECG gating signal seen by the system.
6.	Vertical height indicator 	The Vertical Height Indicator displays the vertical height of the table in relation to ISO center.
7.	Horizontal Cradle Position 	The Horizontal Cradle Position displayed is the position of the cradle based on the established anatomic reference of the patient. This reference is established using the internal or external landmark. The number is preceded by an S if the position is superior to the reference point or an I if the position is inferior to the reference point.
8.	Gantry tilt indicator 	The Gantry Tilt Indicator shows one of two displays. One display, which is the normal display, indicates the current gantry tilt. If the gantry tilt is superior, the number is preceded with an S. If the gantry tilt is inferior, the number is preceded with an I. This is enlarged during Prescribed Tilt button on the Scan Control Interface Module flashes.
9.	Patient Name* 	Patient Name is displayed as input on Patient Information screen. - The system automatically display 2 lines when the name is long. - Font size is automatically changed by the number of characters.
10.	Patient ID* 	Patient ID is displayed as input on Patient Information screen. Font size is automatically changed by the number of characters.
11.	Age, Birthdate and Gender 	Age or Birthdate, and Gender are displayed as input on Patient Information screen. - The patient below one year old, 0 is displayed in Age. - Font size is automatically changed by the number of characters.

SECRET DE AFACERI

No.	Component	Description
12.	Exam Description* Chest-Abdomen	Exam Description is displayed as input on Patient Information screen. Font size is automatically changed by the number of characters.
13.	Default Patient Positioning buttons 	These buttons define the location where the cradle stops with a set pedal. The location can be set by GE Field Engineers.
14.	Tools button 	There are several functions related to the Display setting.
15.	Movie button 	There are some movies that help to explain CT examination procedures to patients.

*Components can Hide/Show during an exam. If you want not to show the Patient Information item during every exam or want to select show components, this can be set by a GE Field Engineer. Contact GE Service to change the selection.

Figure 7-10: A picture of a system with highlighted collision sensor

CONFIDENTIAL

SECRET DE AFACERI

If the system has Emergency Patient on the Gantry display option, the Gantry display includes 2 buttons.

Figure 7-11: Gantry display with Emergency Patient Option and ECG trace on Gantry display option



Table 7-8: Gantry display component with Emergency Patient Option and ECG trace on Gantry display option

No.	Component	Description
1.	ECG trace 	To display ECG trace on Gantry display.
2.	Emergency Patient 	To start an emergency exam.

Figure 7-12: Enabled ECG trace on Gantry display

CONFIDENTIAL

SECRET DE AFACERI

Breath hold lights

Figure 7-19: Breath hold lights


Table 7-13: Breath hold light components

No.	Description
1.	Breathe light
2.	Hold breath light
3.	Elapsed breath hold time in seconds

Internal components

Table 7-14: Gantry internal components

Component	Description
Tube and Col-limator	The anode heat capacity with a Performix™ 40 Plus tube is 7.0 million heat units (MHU) and the cooling rate is 1,070 kilo heat units per minute (KHU/min). Performix™ 40 Plus operates at 48kW or 72kW, depending on configuration.
High Fre-quency Gen-erator	The high frequency generator is composed of an auxiliary box, inverter and high voltage tank. Each tank occupies one cubic foot in the gantry. Together the cathode and the anode provide 72 kilowatts (kW) of the Performix™ 40 Plus tube.
Detector/Data Acquisition System	The detector is comprised of a solid scintillator material known as HiLight. HiLight detectors, provides 98 % absorption efficiency. The V-Res™ Detector enables Volume CT with a true 64-channel design with 40mm coverage.
Slip Ring	Slip Ring provides data communication path from Detector/DAS to the Digital Interface Processor from there to the system disk.

CONFIDENTIAL

SECRET DE AFACERI

Figure 11-23: Real Time Information screen

The screen is titled "Real Time Information". It contains three main sections: "Patient Handling", "Scanning", and "Injecting", each with a large white input field. At the bottom, there are two timing controls: "InterScan Delay" set to "2.5" and "Monitor Scan Time" set to "1.50".

3. View the **display monitor** and track the IV contrast graph.

- Four viewports display from top left to bottom right: baseline image, real time images, graph, time stamp for each real time image.
- A horizontal line indicating the enhancement threshold only appears on the graph if an ROI is placed on the baseline image.
- The tick marks on the SmartPrep graph may not show equal spacing, which is due to time between acquisition rounding issues. The SmartPrep scans are acquired at the correct timing interval, as prescribed.
- When Dynamic Transition is On, the **Dynamic Transition** button is displayed in yellow on Scan Progress screen. If you want to disable Dynamic Transition, click **Dynamic Transition** on Scan Progress screen. Then click **Scan Phase** to start Scan phase manually.
- When Dynamic Transition is On, the time density curve of Transition ROI is displayed in yellow, while the others are in white.

Figure 11-24: Scan Progress screen with Dynamic Transition On

The screen is titled "Scan Progress". It features a timeline at the top with markers from 0 to 60 seconds. Below the timeline, there are four horizontal bars representing different scan phases. At the bottom, there is a row of buttons: "End Exam", "Next Series", "Monitor Phase", "Scan Phase", "Priority Recon", "Image Check", and "Repeat Last Group". On the right side, there is a list of images with columns for "Images" and "Type".

Images	Type
1-1	Real
2-41	Real
1-96	Real
97-192	Real

Scan phase

The Scan Phase is the start of the actual Scan Prescription following the Monitor Phase.

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
SECRET DE AFACERE

AW VolumeShare 7

Improved department productivity.
Increased diagnostic confidence.
Reduced operating costs.
And that's just the beginning...

Today's diagnostic imaging technology provides an amazing breadth of information in a wide variety of ways. However, analyzing it, compiling what's relevant, and reporting your findings can put a serious crimp in your workflow, causing your operating costs to soar, and possibly even affecting your diagnostic confidence. AW VolumeShare 7 offers an elegant solution to your workflow challenges.

Overview

AW VolumeShare 7 is a multi-modality image review, comparison, and processing workstation with simplicity and power at its core. This workstation features 64-bit technology that allows processing of today's increased slice count studies. With significant improvements in workflow processes and clinical features that may help increase your productivity, AW VolumeShare 7 can elevate your practice to a whole new level.

Highlights

- Contemporary, simplified 2D Viewer user interface is harmonized with new Volume Viewer application, new hardware and operating system with expanded memory.
- Supports new applications in CT and Vascular specialties.
- The Postfetch tool retrieves historical exams for comparison.
- Preference Sharing enables exchange of settings and protocols among users.
- The Key Images feature allows you to flag images and create key objects.
- Usage monitoring dashboard tracks and displays exportable analytics per application.
- Supports Chinese, Japanese, Korean and Russian languages for certain applications
- Connects to redundant license servers for high availability

Visit us:
[www.gehealthcare.com/aw/
platforms/aw-workstation/](http://www.gehealthcare.com/aw/platforms/aw-workstation/)

CONFIDENTIAL
SECRET DEAFACER

244

Features

- Simple drag and drop action for networking, media interchange, and filming.
- Multi-tasking capabilities.
- Active directory integration allows enterprise level user authentication. Directories supported include Microsoft® Active Directory® (MSAD) and other LDAP authentication services
- Patient list to manage images from local workstation or remote systems.
- Postfetch feature, designed to optimize your reading workflow by automatically gathering a patient's prior exams according to your pre-set criteria.
- Connectivity with RIS
- Support for Key Images.
- End Review automates routine filming and networking tasks with just one click.
- Search Advantage for fast and easy search of a patient's exam history on PACS or any other DICOM®-compliant device.
- Enhanced Quick Filters of the Patient List filters studies by Modality, Date, End Review status or Exam Description.
- **DICOM CD/DVD/USB Creation Tool.**
- High-performance database management system.
- 2D Viewer for image display, manipulation, annotation, review.
- Integrated Filmer with enhanced flexibility to perform all filming and data exporting tasks.
- Access to a wide variety of applications for greater diagnostic flexibility.

Industry Standards

AW VolumeShare 7 complies with a wide variety of industry standards to facilitate adoption of features and performance improvements as the computing and medical imaging industry evolves.

Computer Industry Standards for Workstation Class Products:

- POSIX (1003.1, 1003.2, 1003.4)
- ANSI C / ANSI Draft Standard C++ / JAVA®

DICOM Conformance Standards:

- DICOM 3.0 Storage Service Class for RT, CT, MR, CR, X-ray (Angio and R&F), Digital X-ray (DX), MG, NM, PET, U/S, Secondary Capture, Secondary Capture Color DICOM Image Objects. (Service Class User (SCU) for image send and Service Class Provider (SCP) for image receive).
- DICOM 3.0 Query/Retrieve Service Class (SCU and SCP).
- DICOM 3.0 Storage Commitment Service Class (SCU).
- DICOM Print (Color and B&W).
- DICOM Media Interchange (CD-R, DVD+R(W)).

Filming Protocols

- DICOM Print (Color and B&W) and Adobe® Postscript (Color and B&W) for supported printers.

Indications for Use

AW VolumeShare 7 is a review workstation, which allows easy selection, review, processing and filming of multi-modality DICOM images from a variety of diagnostic imaging systems. When interpreted by a trained physician, filmed or displayed images on the AW monitor may be used as a basis for diagnosis, except in the case of mammography images.

Regulatory Compliance

- European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/ECUL1950, CSA 950, IEC/EN 60601-1-2:2007, IEC 60950- A12:201
- IEC 60950-1:2005/A1:2009,
- EN 60950-1:2006+A11:2009 +A1:2010+A12:2011,
- UL 60950-1, CSA C22.2 NO 60950-1-07-CAN/CSA CSA C22.2 NO 60950-1-07-CAN/CSA
- EN 980:2008 - Symbols for use in the labeling of medical devices
- EN 1041:2008 - Information supplied by the manufacturer with medical devices
- ISO 14971:2007 / EN ISO 14971:2012 - Medical devices - Application of risk management to medical devices
- IEC / EN 60601-1-6:2010 - Medical electrical equipment
- IEC 62304:2006 / EN 62304:2006 + AC:2008

GE imagination at work

2

CONFIDENTIAL

SECRET DE AFACERI

240

AW VolumeShare 7

The AW VolumeShare 7 workstation gives you powerful capabilities aimed at increasing your diagnostic confidence, reducing your operating costs, and improving your departmental productivity.

Available Configurations:

VolumeShare is a multi-modality advanced visualization workflow solution that helps to enhance diagnostic confidence and productivity. The AW VolumeShare 7 workstation features software that is optimized for 64-bit technology and multicore processor hardware to provide leading edge performance. It includes AW's premier 3D image analysis package, Volume Viewer, with a suite of volumetric visualization and analysis tools for CT, MR, 3D X-ray and PET. In addition, this package includes the multi-modality 2D Viewer with dedicated functions for review of CT, MR, X-Ray, Angio, DX, U/S and PET images. It also includes the Filmer, a multimedia export tool for creating electronic films, filming and exporting to an internal web server, CD, or DVD.

AW VolumeShare 7 is also available for purchase without Volume Viewer for users who require only 2D image applications. Volume Viewer dependent features will not be supported on this configuration.

Modalities

Supported DICOM modalities include:

- Computed tomography (CT),
- Magnetic resonance imaging (MR),
- Radiofluoroscopy (RF),
- Xray angiography (XA),
- Computed radiography (CR),
- Digital radiography (DX),
- Mammography (MG),
- Nuclear medicine (NM),
- Positron emission tomography (PET),
- Ultrasound (US)
- Structured Reporting (SR),
- Key Objects (KO)

Volume Viewer Application

- AW's 3D image analysis package with a broad range of multimodality image analysis tools including multiplanar reformat, curvilinear reformatting, volume rendering and lumen navigation.

- For detailed description of Volume Viewer features, please refer to the Volume Viewer Product Data Sheet.

License Management

- Application Usage Monitor tracks how often a particular license was used, when, and by whom. You can compile and view usage reports for a specific period of time. This helps you understand how your systems contribute to departmental productivity and to plan for future license requirements based on actual usage data.
- For floating licenses, a secondary license server can be assigned for redundancy. If the primary license server becomes unavailable, the secondary can be contacted so that work can continue.

User interface

- Displays on one or two 1280x1024 monitors with scroll or optical mouse.
- Simple drag-and-drop mechanism for networking, media interchange, filming and data export options.
- Automatic help message display briefly describes function currently under the mouse pointer.
- Shortcut keys and programmable window/level function keys for accelerated control.
- Multi-tasking capabilities that allow Patient List, 2D Viewer, Filmer and Volume Viewer to run simultaneously with Fast Switch capabilities between applications.

Security

- Full support of Enterprise Authentication with Microsoft® Active Directory® and Novell® eDirectory™. Users can access their AW Workstation with existing Single Sign On credentials. This may streamline HIPAA compliance.
- If Enterprise Authentication is not used, an administrator-defined UNIX login and password is required per workstation. A lock screen feature enables you to lock the workstation display when leaving it momentarily, preventing unauthorized access to patient data.
- Main actions performed on the AW (save, film, network, export) are flagged and saved to an audit trail which includes username, ID, patient name, exam ID, date, and time.
- Multiple users can have different local UNIX logins, use the same AW workstation, and yet keep their user preferences unique.

- AW VolumeShare 7 is secured by Product Network Filters (PNF), which keeps unauthorized users from accessing open ports over the network. You can configure this feature to provide only specific remote devices access to AW ports.
- While not recommended, AW VolumeShare 7 does support ClamAV® anti-virus protection. This program is not activated by default in order to guarantee optimal performance of the AW workstation. If needed, your GE Healthcare field engineer can activate the program. However, it is the customer's responsibility to provide a safe internet connection from the workstation and to keep virus definitions current.

Workflow management

AW VolumeShare 7 provides features designed to optimize your reading workflow:

- The Patient List provides tools that enable you to sort and filter imaging studies stored on the local workstation or on remote systems.
- Postfetch allows you to retrieve a patient's prior related DICOM exams from a remote DICOM host. Retrieval can be triggered by the arrival of a new patient study on the workstation from the network or supported media.
- Optional integration with RIS systems (Windows® OS only) through command line synchronization interface. This allows your RIS system to automatically select an exam on the AW by passing the patient ID, study UID, or accession number. Together with Softswitch, which allows remote RIS system control from the AW, this feature allows a RIS-driven workflow to be used on the workstation.
- Support for IHE Key Image Notes profile, allowing you to mark certain images in the 2D or 3D viewer as key images. Key images are displayed as a separate series in the exam list and accessed or sent to an IHE KIN-supporting PACS. Specific tags and Key Image Notes may be created using the 2D Viewer application.
- Support for external USB disks to serve as DICOM storage media. The USB disk can be accessed for reading and writing via the media button, the same way you would access a CD or DVD. The maximum number of images that can be stored on a USB device is restricted to 300,000 across all exams. USB media must be formatted using the FAT32 file system. See the AW VolumeShare 7 Basic Display/Viewer/Filmer User Manual for more details on USB storage.
- Preference sharing among users. Shareable preferences include: Filmer layouts, Volume Viewer custom protocols, and viewer preferences.

A standard USB flash drive may also be used to share preferences between workstations.

- One-Touch protocols enable you to define an application or preset protocol to launch automatically based on DICOM elements.
- DICOM Query/Retrieve Storage Class User (SCU) and Storage Class Provider (SCP) provide seamless network integration.
- DICOM Storage Commitment SCU lets you know when exams have been archived successfully on DICOM devices, such as PACS, that support Storage Commitment SCP.

End Review

End Review automates the routine tasks required at the end of reviewing each exam. The "End Review" flag in the patient list allows you to mark exams as "Done" after post-processing has been completed. End Review allows you to automatically perform one or several of these actions:

- Print pages prepared in the Filmer to the default printer (DICOM or Postscript filmers supported), with the option to automatically clear the Filmer after printing.
- Save electronic films in the workstation's attached DICOM database.
- Automatically push the entire exam or only the series created on the AW to one or more remote hosts.

Remote Network Host Management

- AW VolumeShare 7 Patient List provides a control panel to select a remote host or destination for network transfer. Icons signify different device types (acquisition, post processing, PACS, etc.) that are accessible to the workstation via the network.
- You can initiate a network transfer simply by dragging an exam, patient, series, or image(s) to the remote host icon. You can activate a remote browser to display detailed exam information from the remote host, including filtering information if supported.

Search Advantage

- This advanced search engine lets you search a patient's exam history on PACS or other DICOM device with a few mouse clicks. Workflow is improved if the remote device is Direct Connected to the AW so you can launch applications without having to DICOM transfer the exam.

Worklists and Filters

- You can filter the patient list by modality, date, end

GE imagination at work

SECRET INTERNAL

AW CONFIDENTIAL

review status, or exam description. Most filters are available on the remote host patient list as well.

- You can filter the patient list further by choosing one or any of these parameters:
 - Modality.
 - Patient name or patient ID.
 - Exam location (hospital name) and exam description.
 - Series description
 - Date and time of day: today, or specified date or date ranges with a specified time or time ranges.
 - Radiologist's or referring physician's name.
 - Accession number.
- Quick Sort (ascending or descending) and Quick Access (entry field) are available for any of the fields displayed at exam level (patient name, patient ID, exam location and description, radiologist's and referring physician's name, date and time).

Queues

- AW VolumeShare 7 manages three queues: network transfer; media interchange, filming. A menu on the patient list gives you easy access to queue status, and an animated icon gives you continuous network queue status information.

DICOM CD/DVD/USB creation tool

- The DICOM media creation tool offers you the flexibility to save and retrieve from supported CD, DVD, or USB storage devices.
- Add selected exam, series, or images by dragging and dropping to the pre-mastering window.
- CD and DVD usage percent is shown before you launch the writing process.
- You can optionally down sample certain 1024²x512² X-ray angiographic images during CD creation.
- CD/DVD composer lets you manage large amounts of data. If data size is larger than medium capacity, you are prompted to supply additional media storage space.
- You can specify the number of copies created during a CD/DVD save session.
- You can embed a Centricity® DICOM viewer on DICOM CD/DVD so data can be reviewed on a PC running a Windows operating system.
- You get lossless JPEG reading of any DICOM CD/DVD media. Lossless writing is available for X-ray exams.
- CD/DVD drives operate at 16x write speeds and read speeds up to 48x depending on choice of media.
- For limitations on USB storage, refer to the AW VolumeShare 7 Basic Display/Viewer/Filmer User Manual.

Database management

- High-performance database management system classifies data according to patient folder description of the DICOM standard: Patient/Exam; Study/Series;

Sequence/Images.

- Automatically deletes images on first-in first-out (FIFO) basis, which can be toggled on/off.
- Lock Exam feature protects specified exams from deletion.
- Permanent display of available disk space facilitates disk management.
- Study anonymization tool modifies DICOM elements to remove protected health information.

2D Viewer

The 2D Viewer is an application used to display, manipulate, annotate, and review 2D images by a trained physician for diagnostic interpretations. This tool has been designed to provide greater user efficiency, providing a set of intuitive, easy-to-use interface features.

Display customization

Display customization allows the user to manage layouts to display data. Key benefits include:

- Allows modification of number of exams or series displayed to facilitate either single or multi-exam reviews.
- Standard layouts give you flexibility to tailor the image display from 1x1 to 8x8. The current layout is retained if the images are from the same modality and orientation.
- Annotation levels allow selection of the image information fields to display.
- Toolbar customization puts you in control of which buttons are displayed to meet your needs.

Study Navigation

- The navigator lets you assign a series to a view on the fly.
- Cine mode also provides temporal, spatial, or manual playback loops.
- With two exams or series loaded, Cine mode supports a side-by-side display format with synchronized playback loops for more efficient comparisons.

Image Review

- The initial image window and level setting is based on the DICOM header. Once displayed, several methods of adjusting image window and level are provided.
- Provides routine image manipulation features:
 - Flip/Rotate, Zoom, Pan, Magnifying Glass
 - Inverse Video: inverts grayscale color map.
 - Display normal: lets you return the image to its default viewing parameters.

You can access the following features using a single mouse click directly on the image. Having these direct manipulation tools easily at hand

GE imagination at work

gives you fast access with less distraction from your review task by eliminating the need to return to the graphical user interface controls: 2D distance, Angle, Report cursor, Box ROI, elliptical ROI and free-hand ROI. The program continuously updates statistics on the fly.

- The Annotation feature lets you highlight areas of interest by adding text and line/arrow pointers to any image feature.
- The Copy/Paste/Erase feature lets you copy/paste/remove any text or graphic placed on an image.
- Cross-reference indicates the position of the current slice over the scout or localizer or any other non-parallel series.
- Save lets you store a copy of the image as it appears on the screen for future review.
- The Key Images feature allows you to flag images and create key objects.

Batch Filming

The Print Series feature lets you automatically batch film an entire series with a single keystroke.

Advanced X-ray Analysis

- You can apply shutters to X-ray images to focus on specific areas within the image. Shutters are black, opaque overlays with elliptical or rectangular cutouts that you can size as desired. Once a shutter is applied, you can move the underlying image to shift the viewable area. Shutters can also be applied to CT and MR images with the Image Matte feature.
- You can apply several levels of edge enhancement filters to an image.
- **Image Subtraction:** A graphical user interface lets you select a mask and subtraction of all images in an associated sequence. Subtracted images may be saved as a sequence with one mouse click.
- **Landscape:** With this tool, you can introduce a percentage of the mask image into the subtracted image for anatomical reference. A graphical user interface lets you specify the percentage of mask to be reintroduced.
- **Pixel Shift:** Auto pixel shift optimizes a mask to image registration. You can also move a mask manually to optimize subtraction quality in a specific region of interest.
- **Split Pixel Shift:** You can split the screen horizontally or vertically for pixel shift.
- **Maximum/Minimum Opacifications:** You can integrate selected images to provide a resulting Max. Op./Min. Op. image.

Filmer

The integrated Filmer enhances the efficiency of the review station and gives you greater filming and exporting flexibility.

This feature supports two modes: the Mini Filmer mode and Full-Screen mode, which provides the ability for full customization and film layout template creation.

Three key mechanisms give you exceptional flexibility:

1. Free Format Filming
2. DICOM Structured Reporting (SR)
3. Data Export (HTML/PDF and JPEG/PNG/MPEG/AVI/QTVR)

With the Filmer you can easily extract significant images from any AW application (2D Viewer, Volume Viewer, READY View, etc.).

- You can film images individually by dragging and dropping to the on-screen Filmer, or by using the F1 keystroke.
- Multiple Image formatting lets you film multiple images in a single page frame with the F2 keystroke.
- With Film MID you can send multiple images to a single filmer frame with the F3 keystroke.
- Batch Filming is supported from applications that provide that capability (Print Series in the Viewer, Batch Film Protocols in Volume Viewer).

Communication between applications and the Filmer is accomplished in the mini Filmer mode, which provides a minimized footprint. The mini Filmer mode provides the following:

- Compatibility with Batch Filming from Volume Viewer (may be optional in your configuration)
- Store and position images transferred from an application.
- Type of export (film, media, database).
- Rapid switching between Full-Screen Filmer and application.

The resultant electronic films become a quick summary of the patient study and radiological interpretation to be reviewed by clinicians and physicians. They can also be saved as independent files (DICOM SR and Secondary Captures) for teaching purposes.

The Filmer works on the WYSIWYG (What You See Is What You Get) approach, and the electronic film can contain one or several pages with specific layouts for each page or all pages.

A flexible Edit Mode provides the ability to easily add, manipulate, format, or delete images from the film. Images can contain text and graphics from measurements and your annotations, and may be window/leveled, magnified, flipped, rotated, or cine. You can add additional annotation to the image in edit mode.

The Preview Mode displays the film as it will be printed or exported, taking into account the layout applied to each page and the compression level specified for non-DICOM exporting (JPEG/PNG and

CONFIDENTIAL

GE imagination at work

SECRET DE AFACERI
AW Workstation - 4

MPEG/AVI/QTVR).

Printing

AW VolumeShare 7 includes Network DICOM Print (B&W and color). For supported printers, AW VolumeShare 7 includes Network PostScript capability (B&W and color). For PostScript and DICOM printers, any printing format created in the Filmer is supported (e.g. non-square matrix formats for rectangular images such as CT run-offs).

- AW VolumeShare 7 postscript printing has been tested on the following devices: Codonics 1660M, 1660MD or Horizon, Lexmark Optra 1650N, 1855N, SC1275N, C710N, C720N, T612 or T614, Seiko 1720D, Kodak DMI3600, Quantum GL2101HD, declared with film/thick paper or with plain paper, Tally T8106, HP LaserJet and Xerox Phaser
- Digital cameras and analog cameras using a 3M-952 protocol (including DASM interface) are **not** supported on AW VolumeShare 7.

Exporting

Data export is integrated in the Filmer, thus providing all image processing tools needed for multimedia image export. The Filmer exports any electronic film to DICOM SR, PDF/HTML, JPEG, PNG, MPEG, AVI or QTVR format. Non-DICOM data can also be saved on a multi-session CD/DVD. The Data Export capability is intended only for publishing and communication, not for diagnostic purposes. Its simplicity is reflected in the different export mechanisms available:

- CD/DVD removable media
- USB flash drive
- Network HTTP and FTP protocols

Cardiac Review and Export

Processing and review of CT, MR and PET cardiac exams with manual oblique reformatted protocols can be exported as a multi-phase Cine movie that allows referring physicians to review exams in a dynamic mode.

System Components

Workstation Configuration

- HP Z440 Workstation
- HELIOS 6 operating system
- Intel® Xeon® E5-1660 v3 Eight Core 3.0 GHz CPU with 20MB Shared L3 Cache
- 2133 MHz Front Side Bus
- 32GB (4x8GB) DDR4 2133 MHz ECC Registered DIMM
- RAM upgradable to 64GB (8x8GB)
- NVIDIA Quadro NVS 310, 1024MB Graphics card
- **1 x 256GB** Solid State Drive for OS and Applications
- **2 x 512GB** Solid State Drive in RAID-0 for image cache.
- Cache, which is subject to overhead can store approximately:

- 8,300,000 256² uncompressed images OR
- 1,900,000 512² uncompressed images OR
- 475,000 1024² uncompressed images OR
- 100,000 2048 x 2560 uncompressed images
- Internal DVD Writer drive for read/write of DICOM CD/DVD media, read/write of Data Export CD/DVD data and service use (DVD Install)
- Integrated dual Ethernet 10/100/1000 Mbit/s ports.
- **1 USB QWERTY (or regional) Keyboard and Mouse**
- Optional serial port kit for Interventional applications

Footprint

- Height 43.18 cm (17.0 in.)
- Width 16.89 cm (6.65 in.)
- Depth 44.47 cm (17.5 in.)
- Approximate Weight 18 kg (40 lbs)

Operating Environment

- Temperature: +5°C to +35°C
- Humidity: 8% to 85% (relative non-condensing)
- Altitude: 0 to 3100m (10,000 ft.)
- Acoustics: LWAd less than 4.3 Bels
- Shock: 40 G peak, half-sine, 2-3 ms

Non-Operating Environment

- Temperature: -40°C to +60°C
- Humidity: 8% to 90% (relative, non-condensing)
- Altitude: 0 to 9100m (30,000 ft.)

Monitors

- **(2) 19" color Flat Panel LCD monitors**
- DICOM Part 14 factory calibrated
- Native resolution **1280 x 1024** (5:4 aspect ratio)
- Approximate Weight: 7.2 Kg (15.87 lbs)
- AC 100 - 120 V, 200 - 240 V: 50 / 60 Hz

Image Networking

- Standard 10/100/1000 Base-T Ethernet for DICOM
- 1000 Base-T dedicated network for optimal Direct Connect performance
- Protocols supported:
 - DICOM 3.0 Storage SCU/SCP and Query/Retrieve SCU/SCP
 - InSite
 - TCP/IP network layer
 - SdCNet supported to query/retrieve from AWS 3.1 and AW 4.0
- AW VolumeShare 7 does not support the AdvantageNET network protocol.
- AW VolumeShare 7 does not support DICOM images from GE Healthcare Signa™ version 5.x 1.5T MR systems
- AW VolumeShare 7 software is supported on the previous generation AW VolumeShare 5 HP Z800 and Z820 workstations.

CONFIDENTIAL

GE imagination at work

SECRET DE AFACERI
AW Workstation - 5

About GE Healthcare

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services help our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems. Our "healthymagination" vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality around the world. Headquartered in the United Kingdom, GE Healthcare is a unit of General Electric Company (NYSE: GE). Worldwide, GE Healthcare employees are committed to serving healthcare professionals and their patients in more than 100 countries. For more information about GE Healthcare, visit our website at www.gehealthcare.com.

GE Medical Systems SCS
283, rue de la Minière
78530 Buc, FRANCE

©2017 General Electric Company.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. This does not constitute a representation or warranty or documentation regarding the product or service featured. Timing and availability remain at GE's discretion and are subject to change and applicable regulatory approvals. Contact your GE representative for the most current information.

GE, the GE Monogram, Centricity, SIGNA and imagination at work are trademarks of General Electric Company.

Adobe is either a registered trademark or trademark of Adobe Systems Incorporated in the United States and/or other countries.

ClamAV is a registered trademark of Cisco Systems.

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

Intel and Xeon, are either registered trademarks or trademarks of Intel Corporation in the United States and/or other countries.

Microsoft, Active Directory, and Windows, are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Novell and eDirectory are registered trademarks of Novell, Inc. in the United States and other countries.

Oracle and Java are registered trademarks of Oracle and/or its affiliates.

SUSE and SLES are registered trademarks of SUSE LLC in the United States and other countries

All other product names and logos are trademarks or registered trademarks of their respective companies.

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
262

Volume Viewer

GE Healthcare's premium Advanced Visualization and Image Processing platform brings you powerful 3D capabilities to help you visualize and interpret your multimodality imaging data with confidence and ease.

Today's imaging modalities provide a wealth of diagnostic information, but also present a challenge to your Radiology workflow as image volumes continue to increase while your time for reading and reporting decreases.

Volume Viewer provides you with a rich 3D image processing toolset aimed at creating and displaying the views you need with little user input, and streamlining interpretation and reporting by providing the visualization tools you need with minimum clicks.

Now, a refreshed Volume Viewer application with an intuitive, modern user interface helps a wide diversity of users to learn and master the expanding portfolio of sophisticated tools and applications that Volume Viewer provides. This customizable user interface maximizes the real estate allocated for displaying images, and provides access to the users' favorite tools directly from the image viewport. In addition, a customizable toolbar provides one click access to commonly used tools to facilitate a productive review.

Overview

Volume Viewer provides excellent 3D visualization and processing capabilities for reading and comparing CT, MR, 3D X-ray, PET, PET/MR and PET/CT datasets. Volume Viewer also features a broad portfolio of high-performance analysis tools, automating routine tasks and helping to make 3D image processing a stress-free component of your routine workflow.

Volume Viewer is available on VolumeShare 7, a multi-modality advanced visualization workflow solution that helps to enhance diagnostic precision and productivity.

Highlights

- Refreshed user interface with more viewing space for clinical images.
- Simplified and intuitive tools for annotation and measurement.
- Customizable toolbars allowing your favorite tools to be one click away.
- Compare multiple exams from different modalities, at loading and on the fly.
- **3D reformatting, MIP/MPR. High-resolution volume rendering.**
- Easy to use point-and-click tools for segmenting structures of interest.
- Comprehensive, one click or two clicks tracking for any vascular structure (Quick Vessel Trace).
- Interactive Summary Table collecting measurements as they are deposited on the images.
- Adaptive layouts that take advantage of the dual monitor capabilities in landscape and portrait orientations.
- Support of Chinese, Japanese, Korean and Russian languages.¹

Visit us:

www.gehealthcare.com/aw/indications/

CONFIDENTIAL

SECRET DE AFACERI

Features

- **Customizable protocol-driven workflows** for a wide range of clinical applications, designed to help enhance your productivity.
- **Refreshed user interface**, allowing one click access to your favorite tools and more space for displaying images. **Keyboard shortcuts** and **quick access tools** for on-viewport controls.
- **Dual monitor and wide monitor** support allowing more viewing space.
- **Fast automated processing**, with Pre-processing and AutoLaunch capabilities.
- **"Smart Compression"** technology for fast and responsive client performance.
- **"Smart layout"**: Enhanced protocols capability, which adapts the screen layout to the loaded datasets.
- **3D and 4D capabilities**: Real time reformat in oblique planes and interactive multiphase cine of CT, MR, PET.
- **Compare mode** for reviewing multiple exams and modalities. **Dynamic Load** lets you can easily drag and drop new series into the current review session.
- **Multi-modality image fusion** between any 2 different series from same exams.
- **Extensive set of 2D/3D ROI tools** for quantitative measurements and comparisons.
- **Summary Table** extends reporting capabilities.
- **Multiple Export capabilities**, in different formats, including **Key Image Notes** and **STL file**¹.
- **Save state**² allows save and restore to pause and resume the review at a later time.
- **High resolution, real-time rendering modes**: Volume Rendering views; Enhanced rendering of structures with high/low contrast (**creation of MIP/MinIP images**); Navigator view for endoluminal fly-through; lumen and curved views. Creation of **multiple volume rendered models** that may be merged into a single view.
- **Advanced tools**, to take advantage of full 3D capabilities: **AutoSelect**, allowing easy point-and-click segmentation of blood vessels and other structures of interest; **AutoContour** feature for one click contouring of findings and regions of interest on CT, MR and PET images; One or two clicks **Quick Vessel Trace** to analyze all vessels, in curved reformat, lumen, or MPR view.
- **Filters** to enhance and smooth CT images to help improve image quality.
- **Dual energy GE Healthcare CT images** supported with dedicated protocols for review.
- **MR Review** supported with dedicated features: Recognition of PSD name, and specific MR parameters at image loading.
- **SUV measurements** on PET images, for 2D and 3D ROIs, and for iso-contour VOI.
- **Dedicated review protocols for PET/CT and PET/MR images**.
- **Dedicated features for 3D XA images**

System Requirements

- AW Server 3.1 and above and recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048)
- AW VolumeShare7 Workstation and above
- Centricity™ Universal Viewer³

Indications for Use

Volume Viewer is a medical diagnostic software that allows the processing, review, analysis and communication of 3D reconstructed images and their relationship to originally acquired images from CT, MR, X-Ray Angiography and PET Scanning devices. The combination of acquired images, reconstructed images, annotations and measurements performed by the clinician are intended to provide to the referring physician clinically relevant information for diagnosis, surgery and treatment planning.

Regulatory compliance

This product complies with the European Council Directive 93/42/EEC Medical Device Directive as amended by European Council Directive 2007/47/EC.

¹This feature is not available on AW Server 3.1.
²Save State is a Volume Viewer feature allowing the creation of a new series in the database, storing the exam in its current review state with any user-generated post-processing data, for future use.

³AW Server 3.1 is not compatible with Centricity Universal Viewer.

GE imagination at work

© 2016 General Electric Company
All rights reserved. Data subject to change
GE, the GE Monogram, imagination at work, Centricity and AutoBone are trademarks of General Electric company
DOC1460639
GE Medical Systems SCS
283 rue de la Miniere
78530 Buc France

CONFIDENTIAL
SECRET DE AFACERI

Volume Viewer

Detailed Features

- Customizable protocol-driven workflows for a wide range of clinical applications, including support for multiple review steps. You are guided through the review steps by a Navigation bar at the top of the User Interface.

Protocol management and loading with "My Protocols"

- Review protocols created with Volume Viewer can be saved as favorite protocols by clicking "My Protocol". The protocol selection page can be filtered to show "My Protocols" only.
- Standard review protocols are provided and can be launched directly from the AW Server applications menu.
Progressive Load enhances image loading performance. This feature has two modes:
(1) Images are loaded in sequential mode (e.g., Reformat)
(2) Images are loaded in interlaced mode (e.g., 3D/VR)
In either case, image review can begin as soon as the first image is displayed.
- Navigation through series and exams can be performed without exiting to the patient list.


Simple user interface

- The refreshed user interface provides one click access to favorite tools and more space for clinical images.
- The page layout is organized to facilitate navigation through the review steps. Guided protocols consisting of several Review Steps walk the unfamiliar user through a clear workflow of creating and manipulating volumetric data, helping to reduce the learning curve by providing a consistent workflow.
- A customizable toolbar is displayed by default at the top of the user interface, but can be docked at the left, right, or bottom instead, by a simple drag & drop. A user's preferred location is saved for future use.
- The toolbar may be customized with a user's choice of favorite tools for one click access.
- An extensive set of advanced tools are grouped into categories on the toolbar to facilitate quick retrieval when needed.
- Every menu is easily configurable by a simple drag & drop in the Toolbar ... so that they are only one click away for your future reviews!

- To reduce long mouse trips, favorite tools are also accessible from each image via a right mouse menu that may be configured with up to 6 favorite tools.
- Customization of the toolbar and tool behaviors with a user's preferences is provided in a single intuitive Preferences menu on the top navigation bar.

The following mouse modes are available to help you manipulate the images:

- Interactive VR adjustment lets you adjust VR opacity interactively with the mouse as Window Level/Window Width
- Direct Paging – just press the left mouse button
- Free hand rotation on 3D views
- Free hand rotation centered on cursor for oblique views
- Percentage fusion mouse mode to easily change the transparency of objects on fused images. There are several fusion modes available for your use.

- Mag glass mouse mode displays zoomed-up imagery around the cursor.

Dual monitor and wide monitor support

- Volume Viewer provides a comprehensive list of layouts that automatically adapt to the monitor configuration.
- You can use landscape or portraits monitors.

GE imagination at work

Volume Viewer -1
CONFIDENTIAL

SECRET DE AFACERI

Workflow enhancements

Pre-processing⁴ features on AW Server

- Preprocessing is a Volume Viewer feature that automatically performs routine processing tasks and saves the results so they are ready for you. When new exams are transferred to the database, Volume Viewer recognizes user-defined keywords in the Series Description and launches the appropriate pre-processing. The results are stored in a Save State object with the original study so you can load it when you are ready for review.

AutoLaunch⁴ features on AW workstation



- AutoLaunch is a Volume Viewer feature that automatically launches and preloads exams, so they are ready for your review.
- When you are processing a study and a new exam is transferred to the AW Workstation, this feature automatically launches the new exam with an appropriate Volume Viewer protocol in the background. When you are ready, one click in the AutoLaunch window instantly brings up the exam in the Volume Viewer application, eliminating waiting time and extra steps to load the exam into computer memory for processing.
- AutoLaunch is compatible with CT, MR and PET single volume protocols of Volume Viewer.
- When combined with optional applications⁴, this feature gives access to data already preprocessed, that you can review by just clicking on AutoLaunch window.

Smart compression

- The Smart Compression technology automatically displays images with full fidelity when the image is still, then uses the selected compression level for increased interaction speed during user manipulations. This allows for diagnostic reads on full fidelity static images with responsive dynamic display even at low bandwidth. On-image visual indicators notify the user when compression is in effect.

Smart layout

- Volume Viewer General Review protocol presents an adaptive layout for single and multi-modality PACS-like reading of data. This protocol distributes and displays all loaded series of a single exam in an intelligent layout to maximize use of available monitor space. With multiple exams, the protocol automatically detects and distributes "like" series for optimal comparison.

3D and 4D capabilities

- Volume Viewer allows 3D point synchronization in any orientation on the fly.
- Real time reformat in oblique planes is available for your exams, as well as simple and double oblique interactive modes.
- To simplify and improve exam screening, Volume Viewer provides an interactive viewport control that lets you perform the following functions: Continuous or contiguous slice paging; slice thickness adjustment and rendering mode (MIP, MinIP, Average, Volume Rendering).
- Quick access to cross reference display gives you the ability to display the current slice location or all slices, on other planes.
- You can review images in 4D mode, by using Cine tool on all planes and 3D views for CT, MR, and PET multiphase data.

Compare mode

- With Volume Viewer you can load several volumes either from the same exam or from different exams. This is useful for exam comparison (Compare Mode) or for Multiphase examinations (of the liver for example). The Compare Mode is also available for PET/CT examinations and takes advantage of a dual screen configuration.
- Multi-phase CT and MR data are listed as separate phases in the Series Selection screen, allowing you to select all phases or a subset of phase data for review. The Series Selection Panel displays the multiphase data split into arterial and portal phases, for example.
- Dynamic Load, compatible with both Volume Viewer and Fusion protocols, lets you drag & drop 3D volumes from CT, MR, PET, and 3D XA modalities into a desired viewport. Together with the Integrated

GE imagination at work

CONFIDENTIAL Volume Viewer -2

SECRET DE AFACERI

Registration option⁵, Dynamic Load allows new volumes to be registered and loaded on the fly. Save State series can be restored as a separate session. For dual monitor configurations, a Save State series can be displayed on the right monitor together with a current session on the left monitor. For single monitor configurations, a user may simply switch between the two sessions.


- The zoom and pan functions are propagated to all images displayed in the same orientation.

2D/3D ROI tools for quantitative measurements

- Volume Viewer contains a set of basic 2D/3D tools: **distance, angle, report cursor, arrow tool, annotation, free hand ROI tool**, that can be placed on fused images as well.
- You can easily deposit and label measurements.
- All 3D ROIs can be customized and color coded to display statistics computed on different intensity range.
- All measurements are considered bookmarks and you can navigate through the list by using the markers displayed in the slider of each image.

Summary table

- The Summary Table collects measurements and organizes them according to finding and exam date.

- The Summary Table allows you to automatically navigate through findings in the image viewports and can be docked or expanded.

Export capabilities

Volume Viewer contains multiple standard options for exporting the results of a review session:

- Save images to the database, as new series
- Save the significant images as Key Image Note objects in the database. Key Image Notes and End Review allow you to flag images of interest as Key Image Notes (IHE profile) and push them to the archiving system when you exit the application.
- Save State object helps you to save the work and create a new series containing all the post processing data, for future review. It saves the current status of Volume Viewer (3D Model, displays, measurements, annotations, etc), as an additional series in the same exam.
- Save STL file from 3D model, so that you can use it for 3D printing¹.


The following features are available to create the sequence of the images to be exported:

- The Batch allows creating of a sequence of rotating 3D views, or batch of reformatted images.
- The Movie tool creates a comprehensive movie including different rotations, zooms, and pan of the image, which can be exported as DICOM series or mpeg file.
- Quick Export: Exports in a single click a batch of rotations of a 3D View or a full batch of contiguous images at the displayed thickness for 2D images.
- Cardiac Review and Export: Processing and reviewing cardiac exams for CT, MR and PET with manual oblique reformatted protocols can be exported into a multi-phase Cine movie that allows the referring physician to review the exam in a dynamic mode.
- Capture images on the fly on your computer via standard Ctrl+C /Ctrl+V (on AW Server only).

Advanced visualization capabilities

Volume Viewer offers various capabilities to display advanced rendering modes:

- 3D Volume Rendering images.
- **MPR views with different slice thickness. The following rendering modes are available for all Thick Slab: MIP, MinIP, Volume Rendering, Average.**
- Navigator views which display interactive endoluminal views exportable as movies.
- Curvilinear reformatting allowing display of curved, lumen and cross section views for various structures (vessels, spine, etc).

Predefined cut planes

Volume Viewer lets the user define cut planes in order to isolate certain structures in the VR model, in the following 16 configurations:

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GE imagination at work


Multiple VR objects

- Volume Viewer allows you to merge up to 8 Volume Rendered models into a single 3D view.
- You can adjust independently the threshold, colors, and transparency of each VR model.

Auto-center fly through with smart cursor

- In fly-through in air or vessel contrast studies (airways, colon, angiography) you can automatically navigate along the centerline of the structure of interest, simply by pressing the forward or the backward key. You can save the fly-through movie with a simple click.
- In addition to the regular 180° or less Navigator view, Volume Viewer enables a Fish-Eye View with any wide camera angle value from 180° to 360°. This provides a view of structures both in front of and behind the users' virtual location on the same image.

Lumen view

- Lumen View provides an unfolded 3D view of a previously outlined structure. The lumen view can be interactively rotated around the centerline of the structure and set the width and field of view.


Advanced 3D tools

Segmentation tools


- Point and click the AutoSelect segmentation tool to add or remove any continuous structures of interest, including bone and vessel structures for CT and vessel structures for MR. It is also possible to remove the CT couch from images with one click.
- You can also access advanced segmentation tools such as Scalpel, Threshold, Dilate, Erode, Subtraction, Intersection, Addition, Filter Floaters.

AutoContour tool

- Autocontour gives you consistent tools for multi-modality review and comparison: the same contouring tool with the same workflow supports CT, MR, and PET images.
- For CT and MR images, the contouring tool will suggest contour boundaries of structures and let you edit and validate the result.

- For PET images, the contouring tool has been adapted for standard SUV calculations.


GE imagination at work

Quick vessel trace (CT and MR images)


- You can perform vessel analysis with one or two point clicks from any CT or MR Angiography images.
- When you deposit one point, the program automatically extracts the vessel and launches vessel tracking for fast review in curved reformat, cross-sectional, best L-section, lumen, and MRP view.
- You can extend tracking proximally or distally for a full view of the vessel.
- If you detect a lesion, you can quickly analyze it by dropping just one point above and one below the vessel section. The vessel tracking launches automatically.

CT image quality

- Enhancing and smoothing filters are provided to help improve CT image quality for review.
- The lung filter helps enhance contours of images reconstructed in standard mode for excellent visualization of lung structures.

- A smoothing filter reduces noise while maintaining high image contrast:

Dual energy protocols

Three protocols are available to enable review of dual energy images acquired on compatible GE Healthcare scanners:


- **DE Calcium** enhances pixels with ratio $(HU_{80}/HU_{140}) \geq 1.25$
- **DE Uric Acid** enhances pixels with ratio $(HU_{80}/HU_{140}) \leq 1.25$
- **DE Custom** enhances pixels with customized ratio (HU_{Low}/HU_{High})

Review of MR studies

Volume Viewer offers support for review of MR studies, through dedicated protocols that can be selected as My Protocol.

- The MR-specific parameters (PSD name, image weighting, scan plane, etc) are recognized at loading and are used to automatically organize the display by sequence in a specific viewport in the layout.
- Anatomy-based protocols for standardized review (e.g. spine, brain stroke...) thanks to guided layouts with display intelligence to help analyze MR data
- Dedicated MR Cardiac Viewer for single and multiple Cine Review. Enables comparison between cardiac series such as Time Course and Myocardial Delayed Enhancement (MDE) images.
- Supports MR Multi Echo, Multiphase, Diffusion series.
- Supports review of parametric series (e.g., ADC, MTT, and other parametric maps) created in READY View or Functool[®] and enables direct measurement on functional maps.

- 2D and 3D ROI propagation to other phases / series / exams allowing you to easily correlate information from multiple sources.
- Improved MR annotation consistency.
- In combination with Integrated Registration³, enables direct access to MR image registration (inter/intra exam), image fusion and Whole Body MR Review protocols
- Ability to automatically bind several MR axial series corresponding to sequential axial locations into a single series. This may be displayed as sagittal or coronal reformats. This is useful for consolidating multi-stage acquisitions for displaying Whole Body MR scanning.


GE imagination at work

CONFIDENTIAL
SECRET DE AFACERI


SUV measurements on PET images

- Volume Viewer supports Standardized Uptake Values (SUV) for image display and measurement. Several SUV scales are available by right-clicking the unit on-view annotation: SUVbw, SUV lbm, SUVbsa, as well as SUV Peak.
- Window/level presets may be user-defined in SUV.
- When DICOM information necessary for SUV calculation (eg height, weight) is missing from the original images, you can enter it and compute SUV values on the fly.
- The SUV values are available in all the basic 2D/3D ROI tools of Volume Viewer and in the thresholded VOI tool.

PET/CT and PET/MR dedicated protocols

- Volume Viewer contains a list of predefined protocols allowing the review of PET/CT and PET/MR data. The screen layouts contain fused views between the morphological and functional image, 3DMIP rendering of PET images.
- PET 4D protocol allows users to load and display gated and dynamic PET series.
- Factory protocols are customizable for your own review.
- Summing/Reframing tool allows you to Reframe/Rebin series for Dynamic and Gated PET data¹.

Dedicated features for 3D XA images

The optional Innova 3D XR application produces XA images (512x512 images or 256x256 images depending on your settings). This format opens the Volume Viewer world to Innova 3D images:

- Single exam: Innova Navigator, Curved Reformat, Multi Oblique, MPVR 3mm Average, Segment Structure
- Compare exam: Compare Axial, Compare Coronal, Compare Sagittal, Reformat
- Create your layouts and protocols
- The Innova options provide access to Workflow booster tools, easy multi-modality review "on-the-fly" dedicated Innova review protocols with dedicated Review Steps.

Summary of Operation

- Volumetric models are loaded by selecting the exam or series. The user can select a protocol category from an anatomical selector or go directly to a Review Layout. In either case, images are loaded progressively in the background; this gives control to the user in just a few seconds after selecting the images. Selecting a Review Layout launches a volumetric display protocol with predefined layout preferences. Review Layouts may be saved and combined to suit the user's workflow. Selecting a protocol category unlocks a variety of visual protocols that include the layout, threshold, rendering mode and filming formats. Some of these protocols direct the user through the process providing capabilities to interactively view and manipulate the model, increasing productivity and consistency for all modalities.

References

⁴Preprocessing/AutoLaunch is compatible with the following purchasable options: AutoBone™ Xpress, CardIQ Xpress Reveal, CardIQ Xpress Function, Advanced CTC Pro 3D EC, Colon VCAR, CT Perfusion 4D. ColonVCAR is not for sale in the US.

⁵Integrated Registration is a GE Healthcare product, sold separately.

⁶READY View and Functool are GE Healthcare products, sold separately.

GE imagination at work

GE Medical Systems SCS
283 rue de la Minière
78530 Buc France

DOC1460639 rev5

General Electric reserves the right to make changes in specifications and features, or discontinue the product or service described at any time, without notice or obligation. This does not constitute a representation or warranty or documentation regarding the product or service featured. Illustrations are provided for informational purposes, and your configuration may differ.

This information does not constitute legal, financial, coding, or regulatory advice in connection with your use of the product or service. Please consult your professional advisors for any such advice. Operation of GE Healthcare products should neither circumvent nor take precedence over required patient care, including human intervention of healthcare providers. GE Healthcare products and services do not code medical procedures. Accurate coding is the responsibility of the provider or billing professional.

GE, the GE Monogram, imagination at work, Centricity and AutoBone are trademarks of General Electric company.

General Electric Company, by and through its GE Healthcare division.

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281

OncoQuant™

Robust tools for routine Oncology diagnosis, treatment follow-up, reporting and clinical trial management.


Oncology follow-up exams make most of routine studies. Due to the expanding breadth of data associated with today's multi-modality, multiple time-point Oncology studies, reviewing these exams can be time-consuming and labor-intensive. AW's Oncoquant application is designed to help organize and display your Oncology data to facilitate your quick review. Both your routine comparisons and advanced clinical evaluations may benefit from the OncoQuant workflow to navigate through findings and interact with results.

Overview

OncoQuant medical diagnostic software streamlines Oncology reading so you can spend less time retrieving studies and preparing exams, and more time reading and reviewing. A true cross-modality Oncology reading platform, OncoQuant helps you correlate and compare CT, MR, PET/CT, and 3D X-ray data. It automates workflow to facilitate comparisons over time and makes reviewing follow-up exams efficient. The Oncology Review protocol and follow-up wizard simplify your assessment, characterization and measurement of findings on the basis of morphologic criteria.

Highlights

- Automatic multi-modality image registration at loading for two or more exams.¹
- Adaptable workflow supports standard criteria such as RECIST 1.0, 1.1,² WHO³
- Dedicated automatic review protocols helping to identify and load like series.
- Right-click menu workflow measurements at baseline and follow-up.
- Intuitive, interactive Summary Table.
- Single-click Quick Report displaying up to four dates including Baseline, Nadir, Prior, and Current exams.

Visit us:

www.gehealthcare.com/aw/applications/oncoquant/

SECRET DE AFACERI

CONFIDENTIAL

282

Features

- Automatic Registration/Synchronization¹.
- Intelligent review, comparison, and follow up with dedicated protocols.
- Multimodality cases management with no upper limit on the number of exams.
- Advanced Visualization tools:
 - Autocontour tool provides consistent contouring.
 - Lung VCAR algorithm for lung nodules in CT.⁴
 - Capture any finding as a structure of interest using the standard arrow tool.
- Summary table:
 - Baseline and follow up management
 - Drag and drop follow up wizard.
 - Target type selection
 - Customizable Morphological Criterion for routine and research.
- Saves and reloads image and measurements of prior Oncology exams.
- Exports DICOM and non DICOM images of findings, measurements.
- Quick report of the entire review including up to the 4 major dates (Baseline, Nadir, Prior, Current).
- Statistical results export in CSV format for data analysis.

System Requirements

- Available on AW Server 3.1 and above and recommended monitor resolution is up to dual 2MP (1600 x 1200) or a single 3MP (1536 x 2048).
- Available on AW4.7 and above.

Pre-requisite software:

- Integrated Registration is required for Automatic Registration and Synchronization.
- Lung VCAR is required for lung nodule segmentation and Digital Contrast Agent (DCA).

Indications for Use

OncoQuant is medical diagnostic software that allows the processing, review, analysis and communication of 3D reconstructed images and their relationship to originally acquired images from CT, MR, X-Ray Angiography, and PET scanning devices. The combination of acquired images, reconstructed images, annotations and measurements performed by the clinician are intended to provide to the referring physician clinically relevant information for diagnosis, surgery and treatment planning, and follow-up over time.

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References


- ¹ Integrated Registration option required.
- ² RECIST is a guideline to perform clinical trials. Rules and guidelines can be found on the official website <http://www.eortc.org/investigators-preo/recist>
- ³ See WHO criteria in the publication: Measures of Response: RECIST, WHO, and New Alternatives, J Clin Oncol 24:3245-3251.
- ⁴ Lung VCAR option required.

© 2015 General Electric Company.
All rights reserved. Data subject to change.
GE and GE Monogram are trademarks of General Electric Company.
* Trademark of General Electric Company
JB32776XX
GE Medical Systems SCS
283 rue de la Minière
78530 Buc France

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